

Recliners

Florin Recliners

F5017

Mobile Healthcare Recliner Infinite Position Non-Locking Back



Standard Product Information

- > Standard with upholstered arms
- > Standard with LiquiCell technology (*Four gel pads strategically placed in seat*)
- > Seat Foam Density - 2.3 HR (High Resiliency) TB117-2013
- > Back Foam Density - 1.75 HR (High Resiliency) TB117-2013
- > All interior plywood is "California Carb Compliant/Phase II."
- > Easy to remove back to clean for infectious disease control
- > Standard three-4" swivel locking casters and one-4" rigid locking caster
- > Mid-ottoman calf support
- > Infinite position non-locking back (back is not limited to any specific position). Patients can adjust back position to their exact preference through simple pushing on arms and pushing back.

Dimensions

Width	33.5"	Seat Depth (kick-out extended)	40"
Depth	38"	Reclined Depth	76"
Height	47"	Distance from wall	20"
Seat Width	24"	Weight Capacity	500 lbs.
Seat Depth	21"	COM (Yds.)	6
Seat Height	21.5"	Weight	158 lbs.
Back Height	26"	Cubic Feet	26
Arm Height	28"		

Features:

F	Upholstered arms (Standard)	H	Heat (Certified to UL 60601-1)
U	Urethane Arm Caps	HM	Heat & Massage (Certified to UL 60601-1)
W	Wood Arm Caps	S	2-Caster Lock & Steer
Z	Multi Cover		

Options:

CAL	CAL TB133	FT4-1	Large Foot Tray
MBB	Moisture Barrier Back	OTH	Oxygen Tank Holder
MBS	Moisture Barrier Seat	TSL1	Left Side (facing) Tray
AHP	Adjustable Head Rest Pillow	TSR1	Right Side (facing) Tray
4S4	One-4" Swivel Locking Caster	TD3	Dual Side Tray
QR	Quick Release Bracket	LFSTR	Large Right Side (facing) Folding Side Tray
IV	IV Pole and Bracket	SFSTR	Small Right Side (facing) Folding Side Tray
DP	Dual IV Poles	PLFSTR	Padded Large Right Side (facing) Folding Side Tray
FBH	Foley Bag Holder	RSTL	Rotating Side Tray Table (left side facing)
UH	Urine Bag Hook	RSTR	Rotating Side Tray Table (right side facing)
TE	Trendelenburg		

