

# Recliners

## Florin Recliners

F5027

### Mobile Healthcare Rocker Recliner Infinite Position Non-Locking Back



#### Standard Product Information

- > Standard with urethane arm caps
- > Standard with three 3" swivel locking casters and one rigid locking caster
- > Standard with LiquiCell technology (*Four gel pads strategically placed in seat*)
- > Seat Foam Density - 2.3 HR (High Resiliency) TB117-2013
- > Back Foam Density - 1.75 HR (High Resiliency) TB117-2013
- > All interior plywood is "California Carb Compliant/Phase II."
- > Easy to remove back to clean for infectious disease control
- > Mid-ottoman calf support
- > Infinite position non-locking back (back is not limited to any specific position). Patients can adjust back position to their exact preference through simple pushing on arms and pushing back.

#### Dimensions

Width	34.5"	Seat Depth (kick-out extended)	38.5"
Depth	38.5"	Reclined Depth	73.5"
Height	43"	Distance from wall	15.5"
Seat Width	25"	Weight Capacity	375 lbs.
Seat Depth	21"	COM (Yds.)	5.5
Seat Height	19.5"	Weight	160 lbs.
Back Height	23.5"	Cubic Feet	23.5
Arm Height	25"		

#### Features:

U	Urethane Arm Caps (Standard)	H	Heat (Certified to UL 60601-1)
W	Wood Arm Caps	HM	Heat & Massage (Certified to UL 60601-1)
Z	Multi Cover		

#### Options:

CAL	CAL TB133	IV	IV Pole and Bracket
MBB	Moisture Barrier Back	DP	Dual IV Poles
MBS	Moisture Barrier Seat	FBH	Foley Bag Holder
AHP	Adjustable Head Rest Pillow	UH	Urine Bag Hook
3S4	One-3" Swivel Locking Caster	FT5	Regular Foot Tray
QR	Quick Release Bracket	OTH	Oxygen Tank Holder

