


Recliners

Tranquility Recliners

T5025

Healthcare Rocker Recliner Infinite Position Non-Locking Back



Standard Product Information

- > Standard with urethane arm caps
- > Standard with LiquiCell technology (*Four gel pads strategically placed in seat*)
- > Seat Foam Density - 2.3 HR (High Resiliency) TB117-2013
- > Back Foam Density - 1.75 HR (High Resiliency) TB117-2013
- > All interior plywood is "California Carb Compliant/Phase II."
- > Easy to remove back to clean for infectious disease control
- > Mid-ottoman calf support
- > Infinite position non-locking back (back is not limited to any specific position). Patients can adjust back position to their exact preference through simple pushing on arms and pushing back.

Dimensions

Width	36.5"	Seat Depth (kick-out extended)	38.5"
Depth	36.5"	Reclined Depth	74"
Height	43"	Distance from wall	15.5"
Seat Width	25"	Weight Capacity	375 lbs.
Seat Depth	21"	COM (Yds.)	6
Seat Height	18.5"	Weight	139 lbs.
Back Height	24.5"	Cubic Feet	23.5
Arm Height	24"		

Features:

U	Urethane Arm Caps (Standard)	H	Heat (Certified to UL 60601-1)
W	Wood Arm Caps	HM	Heat & Massage (Certified to UL 60601-1)
Z	Multi Cover		

Options:

CAL	CAL TB133	IV	IV Pole and Bracket
MBB	Moisture Barrier Back	DP	Dual IV Poles
MBS	Moisture Barrier Seat	FBH	Foley Bag Holder
AHP	Adjustable Head Rest Pillow	UH	Urine Bag Hook
QR	Quick Release Bracket	OTH	Oxygen Tank Holder

