

Recliners

Tranquility Recliners

T5037

Bariatric Mobile Healthcare Recliner Infinite Position Non-Locking Back



Standard Product Information

- > Standard with urethane arm caps
- > Standard with LiquiCell technology (*Four gel pads strategically placed in seat*)
- > Seat Foam Density - 2.3 HR (High Resiliency) TB117-2013
- > Back Foam Density - 1.75 HR (High Resiliency) TB117-2013
- > All interior plywood is "California Carb Compliant/Phase II."
- > Easy to remove back to clean for infectious disease control
- > Standard three-4" swivel locking casters and one-4" rigid locking caster
- > Mid-ottoman calf support
- > Infinite position non-locking back (back is not limited to any specific position). Patients can adjust back position to their exact preference through simple pushing on arms and pushing back.

Dimensions

| | | | |
|-------------|-------|--------------------------------|----------|
| Width | 40.5" | Seat Depth (kick-out extended) | 40" |
| Depth | 38" | Reclined Depth | 73" |
| Height | 47" | Distance from wall | 19" |
| Seat Width | 29" | Weight Capacity | 750 lbs. |
| Seat Depth | 21" | COM (Yds.) | 7 |
| Seat Height | 21.5" | Weight | 158 lbs. |
| Back Height | 26" | Cubic Feet | 39 |
| Arm Height | 28" | | |

Features:

| | | | |
|---|------------------------------|----|--|
| U | Urethane Arm Caps (Standard) | H | Heat (Certified to UL 60601-1) |
| W | Wood Arm Caps | HM | Heat & Massage (Certified to UL 60601-1) |
| Z | Multi Cover | S | 2-Caster Lock & Steer |

Options:

| | | | |
|-----|------------------------------|--------|--|
| CAL | CAL TB133 | FT4-1 | Large Foot Tray |
| MBB | Moisture Barrier Back | OTH | Oxygen Tank Holder |
| MBS | Moisture Barrier Seat | TSL1 | Left Side (facing) Tray |
| AHP | Adjustable Head Rest Pillow | TSR1 | Right Side (facing) Tray |
| 4S4 | One-4" Swivel Locking Caster | TD3 | Dual Side Tray |
| QR | Quick Release Bracket | LFSTR | Large Right Side (facing) Folding Side Tray |
| IV | IV Pole and Bracket | SFSTR | Small Right Side (facing) Folding Side Tray |
| DP | Dual IV Poles | PLFSTR | Padded Large Right Side (facing) Folding Side Tray |
| FBH | Foley Bag Holder | RSTL | Rotating Side Tray Table (left side facing) |
| UH | Urine Bag Hook | RSTR | Rotating Side Tray Table (right side facing) |

